



# *KineMatch*<sup>™</sup> Custom-Fit Patello-Femoral Replacement




# Indications

**Isolated, end-stage patello-femoral disease**

**Typical Patient:**

- Younger – under 55
- Predominantly female (2:1)



# Kinamed Custom PFR

- Customization provides a perfect fit without bone resection
- No bone cuts, IM rods or sculpting is required with KineMatch
- Four custom cases below show variable anatomy of the trochlea

Four left custom femoral components shown on their CT models



# Kinamed Custom PFR



Lateral aspect



Medial aspect



# Kinamed Custom PFR

## A brief introduction to the KineMatch System

- Custom device based on CT data acquired per Kinamed CT scan protocol
- Provided:
  - CT bone model
  - Custom Femoral Implant
  - Custom Drill Guide
  - Standard 3-peg dome patella
  - Loan of a simple instrument set



# Kinamed Custom PFR

## A brief introduction to the KineMatch System


- Femoral implant has fixation pegs on the backside and is fixed with bone cement (see Surgical Technique Guide).
- Patella implant is an “all poly”, symmetrical, 3-peg design that is also cemented.



# Custom PFR versus TKR

## Custom PFR advantages


- Typically young, active patients – too young for TKR
- Patient functions better on their natural, healthy femoro-tibial articulations.
- PFR is a much less invasive than TKR with 1/3 the morbidity, rehab and recovery time.
- PFR with custom device is a quicker surgical procedure than TKR.
- Look at the X-rays to the left – does it make sense to saw off all the bone required for TKR when the natural femoro-tibial articulations are healthy!?



# Custom PFR versus TKR

## Custom PFR advantages

- Eliminates IM invasion with rod
- Eliminates embolization of fat and marrow caused by rods
- Custom PFR is a bone sparing, temporizing procedure even if the disease progresses to other joint compartments.





# Custom PFR versus Standard PFR

## Custom PFR Advantages



- Custom eliminates bone resection and bone sculpting
- Customization provides a precise fit in the trochlear groove to prevent problems of:
  - patellar catching
  - soft tissue impingement
  - poor patellar tracking and stability.



# Custom PFR versus Standard PFR

## Custom PFR Advantages



- Customization simplifies the surgical procedure.
- Custom drill guide provides precise CT-based fit for creation of peg-holes.



# Custom PFR versus Standard PFR


## Custom PFR Advantages

- Cartilage is removed with a curette.
- Custom femoral component is fixed with bone cement.
- Patella is resurfaced with an all-poly domed implant



# ***KineMatch PFR* Clinical Data**

J. Bone Joint Surg. Am 88:1475-1480, 2006



## **Best clinical results published on any PFR device!**

- 25 PFR in 22 patients
- 16 female, 6 male
- 45 years (23 - 51 years)
- Mean follow-up 73 months

## **Results:**

- 100% Survivorship
- 18 “Excellent” & 7 “Good”
- No revision, loosening, subsequent surgery



## Patello-Femoral Replacement Clinical Results

Author	#	FU (years)	Implant	Result
Ackroyd 2005	306	2 - 5	Stryker	87% not revised and complication-free
Arciero 1988	25	3 - 9	Richards	72% good or excellent, 12% revised
Argenson 1995	66	2 - 10	Medinov	85% not revised
Argenson 2005	66	12 - 20	DePuy	56% not revised
Blazina 1979	57	1 – 3.5	Richards	78% “much improved”
Cartier 1990	72	2 - 12	Richards	85% good or excellent, 8% complications
De Winter 2001	26	1 - 20	Richards	61% good or excellent, 19% reoperations
Kooijman 2003	45	15 - 21	Richards	62% not revised
Krajca 1996	16	2 - 18	Richards	88% good or excellent, 19% reoperations
Lubinus 1979	14	0.5 - 2	Link	“All improved”
Merchant 2004	15	2.2 – 5.5	DePuy	93% good or excellent on ADL scale
Smith 2002	45	0.5 – 7.5	Link	64% good or excellent, 19% revised
Tauro 2001	62	5 - 10	Link	45% “satisfactory”, 28% revised
<b>Sisto 2006</b>	<b>25</b>	<b>2.7 – 9.9</b>	<b>Kinamed</b>	<b>100% good or excellent, No complications</b>