The CE mark is valid only if it is also printed on the product label.

Federal Law (USA) restricts this device to sale by or on the order of a licensed physician.

Devices labeled for single use are intended to be used once only, for a single patient, because they may not perform as intended if they are reused. Reuse may lead to failure of the device to perform as intended.

Do not use if package is damaged.

The presence of this symbol on the product label indicates that the material used for the manufacture of this specific medical device contains Di(2-ethylhexyl) phthalate ("DEHP").

MR SAFETY INFORMATION

Location of Use	MR Safety Information
Within the USA	The CarboJet® System has not been tested for safety and compatibility in the MR environment (it has not been tested for heating, migration, or image artifact in the MR environment).
Outside the USA	The CarboJet® System has been evaluated for safety and compatibility in the MR environment and determined to present a low risk for heating, migration, and image artifact.

INDICATIONS FOR USE
The CarboJet System is indicated for the removal of fluid and particulate debris from bone surfaces for site preparation in orthopedic surgery. Irrigated, sculpted surfaces may be optimized to receive bone cement applied for fixative procedures. CarboJet cleaning is recommended immediately prior to the introduction of bone cement.

PROCEDURES FOR USE: PRIOR TO SURGERY
1. There are several types of outlet fittings used on CO₂ supply cylinders/systems. It is essential to check that the pressure regulator supplied with your CarboJet System is equipped with the corresponding inlet fitting. Kinamed provides the pressure regulator with various types of inlet fittings. Please contact Kinamed Customer Service or your sales representative if you need assistance matching the regulator to the CO₂ source.

2. The CarboJet pressure regulator should be inspected for damage, dirt, dust, oil, or grease. Remove dust or dirt with a clean cloth. Do not attempt to clean the regulator with solvents. If damage is found, return the regulator for repair. As described in Step 1 above, the regulator is used outside of the sterile field and should not be sterilized. The warranty label is affixed to the pressure regulator shows the recommended replacement date. Contact Kinamed Customer Service when the pressure regulator is approaching or reaches the end of its warranty period to discuss replacement options.

3. Secure the CO₂ cylinder to a wall, post, or an appropriate cart so that it can not tip or fall. CAREFULLY open the cylinder valve when the pressure regulator is very small amount for only an instant. Then close the valve quickly to blow out any foreign matter in the valve port.

4. Install the pressure regulator to the CO₂ cylinder/source making sure it is securely fitted. Some regulators come fitted with the appropriate washer. Make sure any required washer is placed before connecting the CO₂ source. Carefully open the cylinder valve and check for audible leaks. When using a CO₂ cylinder tank, check the gauge on the regulator to assess tank pressure. The CarboJet regulators provide measuring functions, as summarized in the table below.

Part Number	Description	Measurement	Degree of precision
25-200-0100	CO ₂ Pressure Regulator with CGA 320 Tank Connector	CO ₂ tank pressure in pounds per square inch (psi) or kilopascals (kPa)	±50 psi or ±345 kPa
25-200-0110	CO ₂ Pressure Regulator with CGA 940 Connector		
25-200-0120	CO ₂ Pressure Regulator with 27mm Connector (Italy)		
25-200-0130	CO ₂ Pressure Regulator with DIN 6 Connector		
25-200-0140	CO ₂ Pressure Regulator with DIN477 FA9 Connector		

Regular delivery pressure to the CarboJet System is factory set to 50 psi (345 kPa). The hospital's Biomedical Engineering department should routinely check the regulator's outlet pressure in accordance with maintenance intervals for similar types of equipment.

5. **VERY IMPORTANT:** After the surgery is completed, be sure to close the cylinder valve completely prior to removal of the regulator from the cylinder.

6. Sterilize the instrument set including the handpiece and all non-disposable nozzles in the CarboJet Sterilization Tray per recommended procedures (see below).

DURING SURGERY

7. Open the pouch containing the sterile tube set and deliver the contents to the sterile field. Using the white quick disconnect fittings at each end of the tube set, connect the female end of the tube set to the handpiece and then pass the male end of the tube set out of the sterile field. Connect the tube set to the CarboJet system. Turn on the CarboJet system and attach the handpiece to the tube set. Turn on the CarboJet system.

8. If sterile disposable nozzles are used, open the pouch containing the sterile nozzle and deliver the contents to the sterile field. Select a nozzle and attach it to the nose of the handpiece via the quick disconnect fitting. Push the nozzle into the fitting until it clicks in place.

9. If using the intramedullary nozzle, first ensure that the stainless steel suction tube is properly assembled on the nozzle and the "Y" tube of the tube set is securely supported on the O-rings at the fitting end of the nozzle. Once the nozzle assembly is complete and the tube set is connected to the handpiece, attach the nozzle to the handpiece.

10. Open the nozzle tip. If this is the first time using the nozzle tip, direct the tip of the nozzle to the prepared bone surface and depress the trigger. Move the nozzle tip away from the sterile field, and then pass the nozzle tip through the corresponding fitting on the handpiece.

11. **VERY IMPORTANT:** With the nozzle tip pointing away from the sterile field, start the gas flow by depressing the trigger on the handpiece for a minimum of five (5) seconds. This step clears the air from the lines and provides a convenient check of all fitting connections.

12. Standard orthopedic surgical procedures are followed for site preparation. Following mechanical shaping and sculpting of the bone bed, saline lavage and suction should be used for initial cleaning of debris and fluids.

13. CarboJet cleaning is recommended as the final step prior to introduction of the bone cement. For most sites, 20 to 30 seconds of CarboJet cleaning is sufficient to remove the debris and liquids. Hold a sponge or towel just ahead of the nozzle tip and gently wash the bone surface with the CarboJet system.

14. **NOTE:** For bone use, the coaxial suction tube must be in place with a suction pump operating at all times for safe and effective cleaning. Long bone cleaning is most effective when working from distal to proximal, as follows: After clearing the tip Step 11 above, insert the nozzle tip into the prepared canal. The nozzle may be rotated back and forth to ensure thorough cleaning.

15. Following CarboJet use, close the valve on the CO₂ tank and purge the CarboJet System by depressing the handpiece trigger. Disconnect the tubing from the regulator and handpiece. Discard the tube set and any disposable nozzles; tube sets and disposable nozzles cannot be re-sterilized and are not re-usable.

MATERIALS

The CarboJet handpiece is made from aluminum alloy with stainless steel components. The tube sets are made from PVC tubing with connectors and components made from vinyl, polycarbonate, Delrin acetate, Buna-N, and stainless steel. The pressure regulator is made from stainless steel and chrome plated brass components. Additional material information is available from the company upon request.

CLEANING & STERILIZATION

All hardware components are fully immersible for cleaning with the exception of the regulator and sterile packaged components. The CarboJet handpiece is made from anodized aluminum alloy. **ENSURE THAT ALL CLEANING AGENTS USED ARE COMPATIBLE WITH THIS MATERIAL.** Re-useable nozzles should be removed from the handpiece for cleaning and the Y tube should be removed from its nozzle. No other disassembly is required or recommended.

16. **NOTE:** For bone use, the coaxial suction tube must be in place with a suction pump operating at all times for safe and effective cleaning. Long bone cleaning is most effective when working from distal to proximal, as follows: After clearing the tip Step 11 above, insert the nozzle tip into the prepared canal. The nozzle may be rotated back and forth to ensure thorough cleaning.

17. The tube sets and nozzles (the status of sterility clearly marked on individual label) have been sterilized by gamma radiation and are sterile unless the pouches have been damaged or opened. The handpiece and re-usable nozzles are provided non-sterile and should be autoclaved per the following validated procedures:

Method	Cycle type	Sterilization temperature	Exposure or Dwell time ^a	Dry time
Steam Autoclave ¹	Pre-vacuum	132°C (273°F) minimum	8 minutes	30 minutes
Autoclave ¹ Double Wrapped ²	Pre-vacuum	134°C (273°F) minimum	6 minutes	30 minutes
Autoclave ¹ Double Wrapped ²	Pre-vacuum	134°C (273°F) minimum	6 minutes	30 minutes

^aThis pre-vacuum sterilization cycle is not considered by the United States Food and Drug Administration (US FDA) to be a standard sterilization cycle. Users should only use sterilizers and accessories (such as sterilization wraps, sterilization pouches, chemical indicators, biological indicators, and sterilization containers) that have been cleared by the US FDA for the selected sterilization cycle specifications (time and temperature).

¹Validated with KimGuard® KC600. (Validated to the following standards: ANSI/AAMI/ISO 17665-1:2006 and ANSI/AAMI ST79-2010.)

²This represents the total exposure/dwell time only and not the total cycle time.

CONTRAINdications

Patients with significant pre-existing cardiopulmonary disorders, including for example a patent foramen ovale (PFO), or who are ASA (American Society of Anesthesiology) Class III or higher, may require careful monitoring by the anesthesiologist during CarboJet use to forestall any unanticipated cardiovascular changes associated with the application of the CO₂.

WARNINGS

1. Only Medical Grade CO₂ gas may be used with the CarboJet device. Use of other gas may result in gas embolism, serious injury, or death.

2. Atmospheric air must be cleared from system prior to use. See instruction #11, above.

3. Only Kinamed CarboJet tube sets may be used to connect the CO₂ source with the CarboJet handpiece.

4. Ensure the cylinder pressure to the tube set does not exceed 50 PSI (345 kPa or 3.45 BAR).

5. The CarboJet pressure regulator contains chemicals including lead, known to be the state of California to cause cancer and birth defects or other reproductive harm.

6. Avoid placing the tip of the CarboJet nozzle into or in close proximity to a venous sinus.

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CAUTION

The CE mark is valid only if it is also printed on the product label.

Devices labeled for single use are intended to be used once only, for a single patient, because they may not perform as intended if they are reused. Reuse may lead to failure of the device to perform as intended.

Do not use if package is damaged.

The presence of this symbol on the product label indicates that the material used for the manufacture of this specific medical device contains Di(2-ethylhexyl) phthalate ("DEHP").

MR SAFETY INFORMATION

Location of Use	MR Safety Information
Within the USA	The CarboJet® System has not been tested for safety and compatibility in the MR environment (it has not been tested for heating, migration, or image artifact in the MR environment).
Outside the USA	The CarboJet® System has been evaluated for safety and compatibility in the MR environment and determined to present a low risk for heating, migration, and image artifact.

INDICATIONS FOR USE

The CarboJet System is indicated for the removal of fluid and particulate debris from bone surfaces for site preparation in orthopedic surgery. Irrigated, sculpted surfaces may be optimized to receive bone cement applied for fixative procedures. CarboJet cleaning is recommended immediately prior to the introduction of bone cement.

PROCEDURES FOR USE: PRIOR TO SURGERY
1. There are several types of outlet fittings used on CO₂ supply cylinders/systems. It is essential to check that the pressure regulator supplied with your CarboJet System is equipped with the corresponding inlet fitting. Kinamed provides the pressure regulator with various types of inlet fittings. Please contact Kinamed Customer Service or your sales representative if you need assistance matching the regulator to the CO₂ source.

2. The CarboJet pressure regulator should be inspected for damage, dirt, dust, oil, or grease. Remove dust or dirt with a clean cloth. Do not attempt to clean the regulator with solvents. If damage is found, return the regulator for repair. As described in Step 1 above, the regulator is used outside of the sterile field and should not be sterilized. The warranty label is affixed to the pressure regulator shows the recommended replacement date. Contact Kinamed Customer Service when the pressure regulator is approaching or reaches the end of its warranty period to discuss replacement options.

3. Secure the CO₂ cylinder to a wall, post, or an appropriate cart so that it can not tip or fall. CAREFULLY open the cylinder valve when the pressure regulator is very small amount for only an instant. Then close the valve quickly to blow out any foreign matter in the valve port.

4. Install the pressure regulator to the CO₂ cylinder/source making sure it is securely fitted. Some regulators come fitted with the appropriate washer. Make sure any required washer is placed before connecting the CO₂ source. Carefully open the cylinder valve and check for audible leaks. When using a CO₂ cylinder tank, check the gauge on the regulator to assess tank pressure. The CarboJet regulators provide measuring functions, as summarized in the table below.

Part Number	Description	Measurement	Degree of precision
25-200-0100	CO ₂ Pressure Regulator with CGA 320 Tank Connector	CO ₂ tank pressure in pounds per square inch (psi) or kilopascals (kPa)	±50 psi or ±345 kPa
25-200-0110	CO ₂ Pressure Regulator with CGA 940 Connector		
25-200-0120	CO ₂ Pressure Regulator with 27		

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