KINAMED[®]

INCORPORATED

www.kinamed.com/eifu

Kinamed Incorporated

820 Flynn Road

Sunnyvale CA, 93012 USA

Tel: 1-800-527-5775, Fax: 1-805-384-2792

www.kinamed.com

EU Authorized Representative:

EC REP

MediMark[®] Europe Srl, 11, rue Emile Zola, BP 2332

F-38033 Grenoble Cedex 2, France

B00040N

Y998-328

©2020 Kinamed[®], Inc.

EN English

Kinamed[®] CarboJet[®] CO₂ Bone Lavage System

Instructions for Use

CAUTION

Rx Only

Federal Law (USA) restricts this device to sale by or on the order of a licensed physician.

MD Medical Device

STERILE R

Tube sets and selected nozzles are sterilized using gamma irradiation.

Do not sterilize.

Devices labeled for single use are intended to be used once only, for a single patient, because they may not perform as intended if they are reused. Reuse may lead to failure of the device to perform as intended.

Do not use if package is damaged.

The presence of this symbol on the product label indicates that the material used for the manufacture of this specific medical device contains Di(2-ethylhexyl) phthalate "DEHP".

MR SAFETY INFORMATION

Location of Use MR Safety Information

Within the USA

The CarboJet[®] System has been evaluated for safety and compatibility in the MR environment and determined to present a low risk for heating, migration, and image artifact.

Outside the USA

The CarboJet[®] System has been evaluated for safety and compatibility in the MR environment and determined to present a low risk for heating, migration, and image artifact.

INDICATIONS FOR USE

Prior To Surgery

1. There are several types of outlet fittings used on CO₂ supply cylinders/systems. It is essential to check that the pressure regulator supplied with your CarboJet System is equipped with the corresponding inlet fitting. Kinamed provides the pressure regulator with various types of inlet fittings. Please contact Kinamed Customer Service or your sales representative if you need assistance matching your CarboJet system to your CO₂ cylinder.

2. The CarboJet pressure regulator should be inspected for damage, dirt, dust, oil, or grease. Remove dust or dirt with a clean cloth. Do not use the regulator if oil, grease, or damage is present. No other maintenance is required for the regulator, except as described in Step 4 below. The regulator is used outside of the sterile field and should not be sterilized. The warranty label that is affixed to the pressure regulators shows the recommended replacement date. Contact Kinamed Customer Service when the pressure regulator is approaching or reaching end of its warranty period to discuss replacement options.

3. Secure the CO₂ cylinder to a wall, post, or an appropriate cart so that it can't tip or fall. CAREFULLY open the cylinder valve a few turns clockwise until the flow begins to form bubbles in the manometer valve.4. If the CarboJet pressure regulator is not connected to the CO₂ cylinder before connecting to the CO₂ source, carefully open the cylinder valve and check for audible leaks. When using a CO₂ cylinder tank, check the gauge on the regulator to assess tank pressure. The CarboJet regulators provide measuring functions, as summarized in the table below.

Part Number Description Measurement Degree of Accuracy

25-200-0100 CO₂-Pressure Regulator with CGA 320 Tank ConnectorCO₂ tank pressure in pounds per square inch (psi) or kilopascals (kPa)

±50 psi or ±345 kPa

25-200-0110 CO₂-Pressure Regulator with CGA 940 ConnectorCO₂-Pressure Regulator with DIN 3002 (Italy)CO₂-Pressure Regulator with DIN 4774 FA ConnectorCO₂-Pressure Regulator with DIN4774 FA Connector

Regulator delivery pressure to the CarboJet System is factory set to 50 psi (345 kPa). The hospital's Biomedical Engineering department should routinely confirm the regulator's outlet pressure in accordance with maintenance intervals for similar types of equipment.

5. **VERY IMPORTANT:** After the surgery is completed, be sure to close the cylinder valve completely prior to removal of the regulator from the cylinder.

6. Steam sterilize the instrument set including the handpiece and all non-disposable nozzles in the CarboJet Sterilization Tray per recommended procedures (see below).

DURING SURGERY

7. Open the pouch containing the sterile tube set and deliver the contents to the sterile field. Using the white quick disconnect fittings at the end of the tube set, connect the female end of the tube set to the handpiece and then pass the male end of the tube set out of the sterile field and connect it to the pressure regulator. Both tubing and fittings are simply pushed into place until an audible click is heard. If the quick disconnect fittings do not "click" into place, repeat the steps to engage the screw lock on the corresponding fitting.

8. If sterile disposable nozzles are used, open the pouch containing the sterile nozzle and deliver the contents to the sterile field. Select a nozzle and attach it to the nose of the handlepiece via the quick disconnect fitting. Push the nozzle into the fitting until it clicks in place.

9. If using the CO₂-Nyloc, first ensure that the stainless steel screw tube is properly assembled on the nozzle and the "Y" end of the tube is securely supported on the O-rings at the fitting end of the nozzle. Once the nozzle assembly is complete and the nozzle is attached to the handpiece, connect a standard sterile suction tube to the "Y" port of the nozzle's steel suction tube.10. Open the valve on the CO₂ cylinder tank. If any noise is heard, check to ensure that the washes between the regulator and cylinder is present and that the regulator fitting is properly tightened on the CO₂ cylinder to avoid leaking.11. **VERY IMPORTANT:** With the nozzle tip pointing away from the sterile field, start gas flow by depressing the trigger on the CO₂ cylinder tank for a minimum of five (5) seconds. This step clears the air from the lines and provides a convenient check of all fitting connections.

12. Standard orthopedic surgical procedures are followed for site preparation. Following mechanical shaping and sculpting of the bone bed, saline lavage and suction should be used for initial cleaning of debris and fluids.

13. CarboJet cleaning is recommended as the final step prior to introduction of the bone cement. For most sites, 20 to 30 seconds of CarboJet lavage should provide adequate cleaning. Direct the tip of the appropriate nozzle at the prepared bone bed, debris and liquids. Hold a sponge or towel just ahead of and directly opposite the tip of the nozzle to collect any debris that will be scattered by the jet. Nozzles may be changed as needed during the procedure.

14. **NOTE:** For long bone, the coaxial tube must be in place with a suction pump operating at all times for an effective cleaning. Long bone cleaning is most effective when working distal to proximal, as follows: After clearing the line per Step 11 above, insert the long bone nozzle into the prepared femur, with the suction pump running, and then depress the CarboJet trigger. With the trigger depressed, slowly draw the nozzle proximally to clear the prepared canal. The nozzle may be rotated back and forth to ensure thorough cleaning.15. Once the nozzle is removed from the CO₂ tank and purge the CarboJet System by depressing the handpiece trigger. Discard the tubing from the regulator and handpiece. Discard the tube set and any disposable nozzles, tube sets and disposable nozzles cannot be re-sterilized and are not reusable.

MATERIALS

The CarboJet handpiece is made from stainless steel components. The tube sets are made from PVC tubing with connectors and components made from vinyl, polycarbonate, Delrin acetate, BuNa-N, and stainless steel. The pressure regulator is made from stainless steel and chrome plated brass components. Additional material information is available from the company upon request.

STERILITY AND HANDLING

The tube sets and selected nozzles (status of sterility clearly marked on individual label) have been sterilized by gamma radiation and are sterile unless the pouches have been damaged or opened.

All instruments and components are provided non-sterile and must be cleaned and sterilized before surgical use. The CarboJet handpiece is made from anodized aluminum alloy. **ENSURE THAT ALL REUSABLE INSTRUMENTS ARE COMPATIBLE WITH THIS MATERIAL.** Re-usable nozzles should be removed from the handpiece for cleaning and the Y-shield should be removed from its nozzle. No other disassembly is required or recommended. Sterilization of instruments must be thoroughly cleaned prior to autoclaving. Inspect devices for dryness prior to sterilization. Autoclave using the following recommended procedures.

Location of Use Method Cycle type Sterilization temperature Exposure or Dwell Time Dry Time

Within the USA Steam Autoclave Double Wrapped¹

132°C (267°F) Minimum

6 Minutes Minimum

30 Minutes Minimum

Outside the USA Steam Autoclave Double Wrapped¹

134°C (273°F) Minimum

6 Minutes Minimum

30 Minutes Minimum

¹Validated with KimGuard[®] KC600. (Validated to the following standards: ISO 17655-1:2006 and ANSI/AAMI ST79-2010.)²This represents the total exposure/well time only and not the total cycle time.

The reusable instruments should be visually inspected for wear and tear, routinely prior to each surgical use. Visual inspection should be performed to look for imperfections or deterioration (such as corrosion, pitting, significant discoloration, cracking, breaking, disassembly, etc.). For reusable instruments that show evidence of such deterioration, the user should first contact Kinamed (or for users located outside the USA, the authorized Kinamed distributor) to determine if repair is possible. If Kinamed determines that repair is possible, the user should thoroughly clean the deteriorated reusable instrument before returning for repair. Kinamed determines that repair is not possible, then the deteriorated reusable instrument should be removed from the handpiece for cleaning and the Y-shield should be removed from its nozzle. No other disassembly is required or recommended. Sterilization of instruments must be accomplished by autoclaving. Instrumentation must be thoroughly cleaned prior to autoclaving.

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