KINAMED
INCORPORATED

www. Kinamed.com/eifu

Kinamed Incorporated
820 Flynn Road, Camarillo, CA 93012 USA
(800) 827-5775 (805) 384-2792 (Fax)
www.Kinamed.com

ISO 13485 Quality Management
certified by BSI under certificate
number FM75124

B00235G
Y998-400

©2024 Kinamed®, Inc.

EU Authorized Representative:

EC

REP

MediMark® Europe Sarl. 11, rue Emile Zola. BP 2332
F-38033 Grenoble Cedex 2. France

EN

English
KINAMED®
SuperCable® Curved Mini Passers
Technique Supplement

CE

2797

The CE mark is valid only if it is also printed on the product label.

This document supplements the SuperCable® Iso-Elastic™ Cerclage System Instructions for Use (document B00109) and Surgical Technique (documents B00110 and B00250).

EXPLANATION OF SYMBOLS	
	Medical device
	Federal Law (USA) restricts this device to sale by or on the order of a licensed physician.
	Single use instruments are sterilized using gamma irradiation.
	Single sterile barrier system with protective packaging inside.
	Do not resterilize
	Do not use if package is damaged and consult instructions for use.
	Do not re-use . Devices labeled for single use are intended to be used once only, for a single patient, because they may not perform as intended if they are reused. Reuse may lead to failure of the device to perform as intended.
	Catalog number
	Lot number
	Caution
	Manufacturer
	Date of Manufacture
	Use-by date
	Unique Device Identifier
	Consult instructions for use or consult electronic instructions for use.

MR SAFETY INFORMATION

Location of Use	MR SAFETY INFORMATION
Within the USA	The SuperCable® implant has not been tested for safety and compatibility in the MR environment (it has not been tested for heating, migration, or image artifact in the MR environment).
Outside the USA	The SuperCable® implant has been evaluated for safety and compatibility in the MR environment and determined to present a low risk for heating, migration, and image artifact.

As stated in the Instructions for Use (B00109), the **INDICATIONS FOR USE** for the SuperCable Iso-Elastic Cerclage System are as follows: (1) Repair of long bone fractures due to trauma or reconstruction; (2) Reattachment of the greater trochanter in total hip arthroplasty, surface replacement arthroplasty, or other procedures involving trochanteric osteotomy; and (3) Sternotomy closure.

As also stated in the Instructions for Use, the **CONTRAINDICATIONS** are as follows: (1) Active or suspected infection, either systemic or localized, in or around the implant site; (2) Patient conditions, mental or neurological, that would tend to impact the patient's ability to follow physician's instructions during the post-operative healing phase; (3) Skeletally immature patients; (4) Compromised vascularity that would inhibit adequate blood supply to the fracture or operative site; (5) Insufficient quality or quantity of bone, which would inhibit rigid device fixation; (6) Any disease affecting the support and fixation of the prosthesis; (7) Distant foci of infection, such as genitourinary, pulmonary, skin or other sites which may spread to the implant site. The foci of infection should be treated prior to, during, and after implantation.

INSTRUCTIONS FOR USE

This Technique Supplement contains information on use of the SuperCable Curved Mini Passers, which are manual surgical instruments that may be useful when using SuperCable under the labeled indications. Refer to Kinamed document B00109 for the SuperCable Iso-Elastic Cerclage System Instructions for Use. General surgical technique instructions are contained in Kinamed documents B00110 and B00250.

PASSER DESCRIPTION

SuperCable Curved Mini Passers are single use manual instruments that facilitate passage of the SuperCable implant ends, such as for passage through rigid or semi-rigid structures like bone, cartilage, or tendon. The Curved Mini Passer shape and cross section are designed to facilitate passage around or through such structures. Curved Mini Passers are designed for use with SuperCable implants only. Refer to the Surgical Technique (documents B00110 and B00250) for descriptions of the SuperCable implant and other associated instruments.

STERILITY

The SuperCable Curved Mini Passers are provided sterile and are intended for single use only. The package should be examined prior to use for possible breaks in the sterile barrier. Refer to the device package label for additional details including material information and expiration date. Refer to the Instructions for Use (document B00109) for sterility information on the SuperCable implant and other associated instruments.

USAGE

Refer to Kinamed document B00110 or B00250 for instructions on use of the SuperCable implant and associated reusable instruments including the tensioning instrument.

To use a SuperCable Curved Mini Passer, remove the instrument from its package and deliver it to the sterile field. Gently insert both free ends of an unused SuperCable implant approximately 1.5 inches (38mm) into the bore of the Curved Mini Passer to create a temporary assembly. Using conventional suture or needle passing techniques, pass the leading end of the Passer / SuperCable assembly through or around the desired structures.

CAUTIONS

- For sternotomy closure, exercise care to avoid injury to or impingement upon the internal mammary artery or other vessels and nerve bundles. Should vessel damage or suspected vessel damage occur, the device should be removed and the vessel repaired.
- The "Sharp" version of the Curved Mini Passer has a sharp tip and should be used for trans-sternal, not peristernal, cable passage. For any Passer with a leading tip, care should be taken to avoid puncture injuries and glove tears. Forceps or needle drivers should be used for holding and passing.

After passage, remove the Curved Mini Passer by gently pulling out both cable ends. The cerclage procedure may then be continued as described in the Surgical Technique (document

B00110 or B00250). For sternotomy closure, the locking mechanism may be positioned at the edge of the intercostal space to minimize the implant profile. The Curved Mini Passer should be discarded at the end of the procedure.

CAUTIONS

- Do not implant a Curved Mini Passer.
- Avoid wrapping the cables over sharp metal or bone graft surfaces.
- The clasp should be placed in a region of bone that maximizes the conformity between the clasp and underlying surface (bone or allograft). The angle the cable makes with the clasp as it exits the clasp should be as small as is practicable.
- Cable tension should be equalized to the extent possible. When the two free cable ends are inserted into the tensioning instrument, the ends should be pulled taut so as to equalize their length.
- Recommended tensioner settings are meant to assist the surgeon in optimizing performance of the system, not to replace the surgeon's judgment. Care should be taken to control tension in patients with poor bone quality and ideal tension may vary with bone quality or geometry. Reduced bone quality may warrant a lower tension.
- Do not tension cable such that the line on the knob passes the second solid line (H1), exceeding 120 lbs. (530N) of compressive force. Grasp only the knurled portion of the knob and slowly turn while reading the tension level. The indicator marks should be read while turning the outer knob.

If this device is used in the European Union and a serious incident occurs during the use of this device or as a result of use of this device, the incident should be reported to Kinamed and/or its authorized representative and to your national authority.

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