ISO 13485 FM 75124

820 Flynn Road, Camarillo CA, 93012 USA
Tel: 805-384-2748, www.kinamed.com

© 2010

CarboJet® CO₂ Bone Lavage System Instructions for Use

EU Authorized Representative: MediMark Europe
11, rue Emile Zola - B.P. 2332, F-38033 Grenoble Cedex 2 France
Tel +33(0)476864322 Fax +33(0)4761719 82

B00040 Rev. G

CAUTION

Federal Law (USA) restricts this device to sale by or on the order of a licensed physician.

Devices labeled for single use are intended to be used once only, for a single patient, because they may not perform as intended if they are reused. Reuse may lead to failure of the device to perform as intended.

The presence of this symbol on the product label indicates that the material used for the manufacture of this specific medical device contains Di(2-ethylhexyl) phthalate "DEHP".

INDICATIONS FOR USE

The CarboJet System is indicated for the removal of fluid and particulate debris from bone surfaces for site preparation in orthopedic surgery. Irrigated, sculpted surfaces may be optimized to receive bone cement applied for fixative procedures. CarboJet cleaning is recommended immediately prior to the introduction of bone cement.

The CarboJet tube set and disposable nozzles are indicated for single use only.

PROCEDURES FOR USE: PRIOR TO SURGERY

1. There are several types of outlet fittings used on CO₂ supply cylinders/systems. It is essential to check that the pressure regulator supplied with your CarboJet System is equipped with the corresponding inlet fitting. Kinamed provides the pressure regulator with various types of inlet fittings. Please contact Kinamed Customer Service or your sales representative if you need assistance matching the regulator to the CO₂ source.

2. The CarboJet pressure regulator should be inspected for damage, dirt, dust, oil, or grease. Remove dust or dirt with a clean cloth. Do not use the regulator if oil, grease, or damage is present. No other maintenance is required. The regulator is used outside of the sterile field and should not be sterilized.

3. Secure the CO₂ cylinder to a wall or post, or in an appropriate cart so that it can not tip or fall. CAREFULLY open the cylinder valve a very small amount and for only an instant and then close quickly to blow out any foreign matter in the valve port.

4. Install the pressure regulator to the CO₂ cylinder/source making sure it is securely fitted. Carefully open the cylinder valve and check for audible leaks. Check the gauge on the regulator to assess tank pressure. The gauge reads tank pressure only. Regulator delivery pressure to the CarboJet System is factory set to 50 psi (345 kPa).

5. VERY IMPORTANT: After the surgery is completed be sure to close the cylinder valve completely prior to removal of the regulator from the cylinder.

6. Steam sterilize the instrument set including the handpiece and all non-disposable nozzles in the CarboJet Sterilization Tray per recommended procedures (see below).

DURING SURGERY

7. Open the pouch containing the sterile tube set and deliver the contents to the sterile field. Connect the tube set to the handpiece and to the pressure regulator using the white quick disconnect fittings at each end. The tubing end with the female fitting attaches to the handpiece while the end with the male fitting attaches to the regulator. Both tubing end fittings are simply pushed into place until an audible click is heard.

8. If sterile disposable nozzles are used, open the pouch containing the sterile nozzle and deliver the contents to the sterile field. Select a nozzle and attach it to the nose of the handpiece via the quick disconnect fitting. Push the nozzle into the fitting until it clicks in place.

9. If using the intramedullary nozzle, first ensure that the stainless steel suction tube is properly assembled on the nozzle and the "Y" end of the tube is securely supported on the O-rings at the fitting end of the nozzle. Once the nozzle assembly is complete and the nozzle is attached to the handpiece, connect a standard sterile suction tube to the "Y" port of the nozzle's steel suction tube. Connect the opposite end of the suction tubing to the suction tip in the sterile field, connected in turn to a suction canister and pump.

10. Open the valve on the CO₂ tank. If any hissing is heard, check to ensure that the regulator fitting is properly tightened on the CO₂ cylinder to avoid leaking.

11. **NOTE/VERY IMPORTANT:** With the nozzle tip pointing away from the sterile field, start the gas flow by depressing the trigger on the handpiece for a minimum of five (5) seconds. This step clears the air from the lines and provides a convenient check of all fitting connections.

12. Standard orthopedic surgical procedures are followed for preparation of the implant site. Following mechanical shaping and sculpting of the bone bed, saline lavage and suction should be used for initial clearing of debris and fluids.

13. CarboJet cleaning is recommended as the final step prior to introduction of the bone cement. For most sites, 20 to 30 seconds of CarboJet lavage should provide adequate cleaning. Direct the tip of the appropriate nozzle at the prepared bone surface and depress the trigger. Move the nozzle tip as needed to clean the bone of fat, debris, and liquids. Hold a sponge or towel just ahead of the nozzle's tip to collect debris that will be scattered by the spray, particularly when using the wide-angle nozzle. Nozzles may be changed as needed during surgery without disconnecting the CO₂ line.

14. **NOTE:** For long bone use, the coaxial suction tube must be in place with a suction pump operating at all times for effective cleaning. Long bone cleaning is most effective when working from distal to proximal, as follows: After clearing the line per Step 11 above, insert the long bone nozzle into the prepared femur, with the suction pump running, and then depress the CarboJet trigger. With the trigger depressed, slowly draw the nozzle proximally to clean the prepared canal. The nozzle may be rotated back and forth to ensure thorough cleaning.

15. Following CarboJet use, close the valve on the CO₂ tank and purge the CarboJet System by depressing the handpiece trigger. Disconnect the tubing from the regulator and handpiece. Discard the used tube set and any disposable nozzles; tube sets and disposable nozzles cannot be re-sterilized and are not re-useable.

MATERIALS

The CarboJet handpiece is made from aluminum alloy with stainless steel components. The tube sets are made from PVC tubing with connectors and components made from vinyl, polycarbonate, Delrin acetal, Buna-N, and stainless steel. Additional material information is available from the company upon request.

CLEANING & STERILIZATION

All hardware components (unless provided sterile packaged) are fully immersible for cleaning. The CarboJet handpiece is made from anodized aluminum alloy: **INSURE THAT ALL CLEANING AGENTS USED ARE COMPATIBLE WITH THIS MATERIAL.** Re-usable nozzles should be removed from the handpiece for cleaning and the Y tube should be removed from its nozzle. No other disassembly is required or recommended.

Soak the components in Klenzyme™ or equivalent for two minutes, and rinse thoroughly under running tap water. Clean with Manuklenz™ or equivalent using a soft bristle brush. Rinse thoroughly under running tap water.

The tube sets and selected nozzles (status of sterility clearly marked on individual label) have been sterilized by gamma radiation and are sterile unless the pouches have been damaged or opened. The handpiece and re-usable nozzles are provided non-sterile and should be autoclaved per the following validated procedures:

Prevacuum: Wrapped, 4 minutes, 132 °C minimum

Gravity: Wrapped, 30 minutes, 15 psi (103 kPa), 121 °C minimum

Gravity: Unwrapped, 15 minutes, 15 psi (103 kPa), 121 °C minimum

The following cycles are also validated per EN554 for use in the European Community:

Prevacuum: Wrapped, 121° C for 15 minutes

Prevacuum: Wrapped, 134° C to 137° for 3 to 3.5 minutes

CONTRAINDICATIONS

Total joint arthroplasty patients with significant pre-existing cardiopulmonary disorders may require careful monitoring by the anesthesiologist during CarboJet use to forestall any unanticipated cardiovascular changes associated with the application of the CO₂.

WARNINGS

1. Only CO₂ gas may be used with the CarboJet device. Use of other gas may result in gas embolism, serious injury, or death.

2. Atmospheric air must be cleared from system prior to use. See instruction #11, above.

3. Only Kinamed CarboJet tube sets may be used to connect the CO₂ source with the CarboJet handpiece.

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Sistema de lavado do hueso com CO2 CarboJet®
Instrucciones de uso
Representante europeo autorizado: MediMark Europe
11, rue Emile Zola - B.P. 2332, F-38033 Grenoble Cedex 2, Francia
Tel +33 (0) 4 76 86 43 22 Fax +33 (0) 4 76 17 19 82
CAUTION
 Federal Law (USA) restricts this device to sale by or on the order of a licensed physician.
 Devices labeled for single use are intended to be used once only, for a single patient, because they may not perform as intended if they are reused. Reuse may lead to failure of the device to perform as intended.
 The presence of this symbol on the product label indicates that the

material used for the manufacture of this specific medical device contains Di(2-ethylhexyl) phthalate "DEHP".

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CarboJet® CO2 骨洗浄システム
使用説明書
ヨーロッパ認可代理店：MediMark Europe
11, rue Emile Zola - B.P. 2332, F-38033 Grenoble Cedex 2, France
電話 +33 (0) 4 76 86 43 22 ファックス+33 (0) 4 76 17 19 82
CAUTION
 Federal Law (USA) restricts this device to sale by or on the order of a licensed physician.
 Devices labeled for single use are intended to be used once only, for a single patient, because they may not perform as intended if they are reused. Reuse may lead to failure of the device to perform as intended.
 The presence of this symbol on the product label indicates that the

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CarboJet® CO2 气体洁骨系统
使用说明书
欧洲授权销售代理：MediMark Europe
11, rue Emile Zola – B.P. 2332, F-38033 Grenoble Cedex 2, 法国
电话：+33(0)4 76 86 43 22， 传真：+33(0)4 78 17 19 82
CAUTION
 Federal Law (USA) restricts this device to sale by or on the order of a licensed physician.
 Devices labeled for single use are intended to be used once only, for a single patient, because they may not perform as intended if they are reused. Reuse may lead to failure of the device to perform as intended.
 The presence of this symbol on the product label indicates that the

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