

CarboJet CO₂ Bone Preparation System

25-100-2001 Carbojet CO₂ Tubeset, Sterile Packed (each)



Standard Instruments

25-200-0200 CarboJet Handpiece



25-200-0220 CarboJet Angled Tip Nozzle (peg holes, shoulder)



25-200-0230 CarboJet Wide-Angle Nozzle (knee)



25-200-0242 CarboJet Femoral Canal Suction Tube (12mm dia.)




25-200-0244 CarboJet Femoral Canal CO₂ Nozzle




25-200-0246 CarboJet 40 Degree Nozzle (knee, glenoid)



25-200-0300 CarboJet Sterilization Tray



25-200-0110 CarboJet Pressure Regulator



(with CGA 940 pin-index yoke tank connection)

Surgeon Testimonials

"I use it on Uni Knees. We often get a lot of bone bleeding – even with a tourniquet during cementation. The CarboJet dries the bone superbly for implantation. Outstanding addition to my surgical technique!!"

Airell Nygaard, MD, Yosemite Joint Replacement, Sonora, CA, USA

"CarboJet provides my patients with a much cleaner, drier bone bed allowing for better cement penetration and hence a better mechanical interface between bone and implant for secure long-term fixation."

Richard "Dickey" Jones, MD, Dallas, TX, USA

"I have made gas jet lavage with CarboJet the critical last step in bone preparation in all my cemented arthroplasty cases. The removal of additional marrow elements that could otherwise form embolic debris during cement pressurization is important to patient safety."

H.M. "Mac" Reynolds, MD, Oakland, CA, USA, In Memoriam



KINAMED
INCORPORATED

Quality Care. Clinically Proven.

For more information:

Phone (805) 384-2748

Toll-Free (800) 827-5775

Fax (805) 384-2792

Website www.kinamed.com



820 Flynn Road, Camarillo, CA 93012-8701



CarboJet®

CO₂ Bone Preparation System



A Better Way to Clean Bone

A **clinically proven** system that removes lipids/marrow elements and fluids from the cement interface to improve cement penetration and bond strength.

Increase Cement Penetration^{1,2}


Increase Bone-Cement Interface Strength³

Reduce Opportunity for Micro-Emboli⁴

Essential for Tourniquetless TKA^{5,11}

Reduce OR Time and Cost⁶



Tibia After Resection




After Pulsatile Saline Lavage



After CarboJet®



Knee and Uni



Hip



Shoulder

Nozzles are available for use in TKA, UKA, THA, TSA, OCA and other cemented reconstructive applications.



KINAMED
INCORPORATED

Quality Care. Clinically Proven.



The **Gold** Standard for Cement Bed Preparation

Concerned About Loosening?

Aseptic loosening is now recognized as the predominant failure mode in primary knee arthroplasty.^{7,8} A recent study⁷ of 938 primary knees found that aseptic loosening was the most common mechanism of failure, accounting for 28% of all failures. In UKA, aseptic loosening accounted for more than 60% of failures.⁷ These data underscore the importance of achieving optimal cement fixation via meticulous cleaning and drying of the bone bed.

A Simple Solution

CarboJet's CO₂ gas jet quickly and thoroughly cleans and dries the bone bed by bringing blood, saline and, most importantly, lipids/fatty marrow elements to the surface where they are easily collected and removed. Cleaning and drying with CarboJet takes no more time than is typically required for drying with lap sponges. The superior cleaning achieved with CarboJet improves cement penetration depth^{1,2} thus offering the potential for a reduction in bone-cement interface stress⁹ and increased cement mantle toughness.¹⁰ Increased bone-cement interface strength is the result.³



Saline lavage removes gross debris, but fluid remaining in the interstices of bone prevents thorough cleaning.




A compressed gas jet displaces fluid and fluid-suspended debris, thereby cleaning and drying porous structures.

Clinically Proven

CarboJet has been proven to be safe and effective in multiple clinical studies and in tens of thousands of joint reconstructions. Discover why so many surgeons are making CarboJet a standard part of their cement technique.

Cement mantle toughness correlates with cement penetration depth.¹⁰



Bone-cement interface strength is 58% higher with CarboJet cleaning versus saline lavage.³




1. Goldstein (2007) **Improvement of cement mantle thickness with pressurized carbon dioxide lavage.** ISTA. Paris, France.
2. Woodgate (2008) **A radiological comparison of cement mantle thickness around TKA.....**Australian Orthopaedic Association.
3. Ravenscroft, Stanley et al (2010) **Bone-Cement Interface Strength in Distal Radii Using Two.....**Hand Surgery 15(2):95-98.
4. Lassiter, Bolognesi et al (2010) **Intraoperative embolic events during TKA with use of pulsatile.....**ORS. New Orleans, USA.
5. Jones (2011) **Total Knee Arthroplasty without the use of a tourniquet.** Seminars in Arthroplasty 22:176-178.
6. Stiehl (2014) **Mechanical Performance of a Self-Unplugging Surgical Suction.....**Reconstructive Review 4(1):18-22.
7. Berend, Lombardi, Barnes, Bolognesi, Ritter et al (2013) **Mechanism of Primary Knee Arthroplasty.....**AAOS.
8. Lewold et al (1998) **Revision of unicompartmental knee arthroplasty: outcome in.....**Acta Orthop Scand.69(5):469-74.
9. Thompson et al (2010) **The Importance of a Good Cement Mantle with an All-Poly Inlay UKA.** ORS. New Orleans, USA.
10. Graham et al (2003) **Effect of Bone Porosity on the Mechanical Integrity.....**J Bone Joint Surg Am. 85:1901-1908.
11. Meneghini (2018) **Tourniquetless TKA...Decreases Pain and Opioid Consumption in Females.....**AAOS. New Orleans, USA.

Improved Cement Mantle



CarboJet enables cleaning and drying of posterior aspects of tibial and femoral surfaces.



Bone bed prepared with pulsatile saline lavage. Arrows indicate radiolucent lines.




Bone bed prepared with syringe saline lavage and CarboJet CO₂ Bone Preparation System.

Essential for Tourniquetless TKA

Tourniquetless Total Knee Arthroplasty with Modern Perioperative Protocols Decreases Pain and Opioid Consumption in Females

A recent study,¹¹ presented at AAOS 2018 by R. Michael Meneghini, M.D., demonstrates reduced opioid consumption in female patients when a tourniquet was not used in total knee arthroplasty (TKA). CarboJet was used to maximize cement interdigitation – an essential component of the tourniquetless TKA technique. The ability to reduce opioid consumption using this technique is significant, since opioid addiction has become a national crisis, and protocols that result in reduced patient pain are being sought.



Scan for Video of CarboJet in Tourniquetless TKA:



B/C FPO.

Do NOT print red dieline.