Product Insert
Hoja explicativa del producto
Folheto do produto
Notice explicative
Nota esplicativa
Packungsbeilage

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NeuroPro® B00016 G

KINAMED® NeuroPro® System

MATERIALS

1. Screws: Ti 6Al 4V-ELI, ISO 5832-3, ASTM F-136
2. Plates and Panels: Ti CP, ISO 5832-2, ASTM F-67

INDICATIONS

1. Internal fixation of fractures and osteotomies of the craniofacial skeleton.
2. Internal fixation of cranial bone flap osteotomies.
3. Reconstruction of bony defects and deficits in the craniofacial skeleton.

CONTRAINDICATIONS

1. Active or suspected infection, either systemic or localized, in or around the implant site.
2. Patient conditions, mental or neurological, that would tend to impact the patient's ability to follow physician's instructions during the post-operative healing phase.
3. The Kinamed NeuroPro® System is not indicated for use in the spine or high load bearing applications.
4. Demonstrated sensitivity to Titanium or its alloy (Ti6Al4V-ELI).

WARNINGS

1. An implant should never be re-used. Even though an implant may appear undamaged, previous handling and in-service stresses may have created imperfections that would reduce the service life of the implant.
2. Repeated bending back and forth of the plates should be avoided as this action may weaken the plate leading to plate breakage.
3. In the absence of functional, healed bone, implant failure may occur if the bone segments are subjected to repetitive loading over time.
4. Patients should be questioned regarding sensitivity to metals. If there is a question pertaining to the patients tolerance for titanium or its alloys appropriate testing should be performed.
5. Pilot hole drill bits should be used for a single surgical case only.

PRECAUTIONS

1. Once removed from their original packaging, implants and instruments should be stored and autoclaved in the organizer trays provided in the system to prevent contact with items of dissimilar metals.
2. Adequate inventory of the various sizes and configurations of implants should be available in the organizer tray at the time of surgery to meet the requirements of each specific surgical case.
3. Following use instruments should be thoroughly cleaned prior to being replaced in the organizer tray for sterilization.
4. Drills should be inserted only into the instruments for which they are labeled.

STERILITY AND HANDLING

All items in the system, including the implants, are supplied **non-sterile**. Sterilization of the instruments and implants is accomplished by autoclaving. The following are recommended procedures:

Note: Sterilization times given below represent exposure time only and not total cycle time.

Pre-Vac Autoclave
132° C, +/- 1° C (270° F, +/- 2° F)
15 minutes wrapped

Caution: Federal Law (USA) restricts this device to sale by or on the order of a licensed physician.

European Authorized Representative:
MediMark Europe, Zola B.P. 2332, 38033 Grenoble Cedex 2, France

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© 2007