www.kinamed.com/symbols_glossary

Kinamed Incorporated
820 Flynn Road, Camarillo, CA 93012 USA
(800) 827-5775 (805) 384-2792 (Fax)
www.kinamed.com

B00016K

©2019 Kinamed®, Inc.


EU Authorized Representative:

MediMark® Europe Sarl. 11, rue Emile Zola. BP 2332 F-38033 Grenoble Cedex 2. France

English
KINAMED® NeuroPro® System
Instructions for Use

CAUTION

Federal Law (USA) restricts this device to sale by or on the order of a licensed physician.

Devices labeled for single use are intended to be used once only, for a single patient, because they may not perform as intended if they are reused. Reuse may lead to failure of the device to perform as intended.

Do not use if package is damaged.

MR SAFETY INFORMATION

Location of Use	MR Safety Information
Within the USA	The NeuroPro® implant has not been tested for safety and compatibility in the MR environment (it has not been tested for heating, migration, or image artifact in the MR environment).
Outside the USA	The NeuroPro® implant has been evaluated for safety and compatibility in the MR environment and determined to present a low risk for heating, migration, and image artifact.

MATERIALS

- Screws: Ti 6Al 4V-ELI, ISO 5832-3, ASTM F-136
- Plates and Panels: Ti CP, ISO 5832-2, ASTM F-67

INDICATIONS

- Internal fixation of fractures and osteotomies of the craniofacial skeleton.
- Internal fixation of cranial bone flap osteotomies.
- Reconstruction of bony defects and deficits in the craniofacial skeleton.

CONTRAINDICATIONS

- Active or suspected infection, either systemic or localized, in or around the implant site.
- Patient conditions, mental or neurological, that would tend to impact the patient's ability to follow physician's instructions during the post-operative healing phase.
- The Kinamed NeuroPro® System is not indicated for use in the spine or high load bearing applications.
- Demonstrated sensitivity to Titanium or its alloy (Ti6Al4V-ELI).

WARNINGS

- An implant should never be re-used. Even though an implant may appear undamaged, previous handling and insertive stresses may have created imperfections that would reduce the service life of the implant.
- Repeated bending back and forth of the plates should be avoided as this action may weaken the plate leading to plate breakage.
- In the absence of functional, healed bone, implant failure may occur if the bone segments are subjected to repetitive loading over time.
- If extremely dense bone is encountered while using Quick Tap self-drilling screws, it is recommended that a pilot hole be drilled using a drill bit from the system to prevent possible screw breakage during insertion.
- Drill bits provided in the system must be used to assure correct pilot hole size.
- Patients should be questioned regarding sensitivity to metals. If there is a question pertaining to the patient's tolerance for titanium or its alloy appropriate testing should be performed.
- Pilot hole drill bits should be used for a single surgical case only.
- Once screw-head is fully seated into plate, avoid over-tightening to prevent possible screw breakage or pressure necrosis of the bone under the plate.

PRECAUTIONS

- The 4mm screws require less torque to seat fully in dense bone due to less frictional resistance over their shorter length.
- The 5mm screws are used in soft bone where more resistance to stripping is required and the bone thickness is great enough to prevent penetration of the screw through the inner cortex or "table" of the skull.
- Once removed from their original packaging, implants and instruments should be stored and autoclaved in the organizer trays provided in the system to prevent contact with items of dissimilar metals.
- Adequate inventory of the various sizes and configurations of implants should be available in the organizer tray at the time of surgery to meet the requirements of each specific surgical case.
- Following use, instruments should be thoroughly cleaned prior to being placed in the organizer tray for sterilization.
- Drills should be inserted only into the instruments for which they are labeled.

STERILITY AND HANDLING

All items in the system, including the implants, are supplied **non-sterile** and must be sterilized prior to surgical use. Sterilization is accomplished by autoclaving. Instruments must be thoroughly cleaned prior to autoclaving. Inspect devices for dryness prior to sterilization. Autoclave using the following recommended procedures:

Note: Sterilization times given below represent exposure time only and not total cycle time.

Location of Use	Method	Cycle Type	Sterilization Temperature	Full Cycle Time	Dry Time
Within the USA	Steam Autoclave (Double wrapped in 1-ply polypropylene wrap.)	Pre-Vacuum	132° C (270° F) Minimum	15 Minutes Minimum	80 Minutes Minimum
Outside the USA	Steam Autoclave (Double wrapped in 1-ply polypropylene wrap.)	Pre-Vacuum	134° C (273° F) Minimum	12 Minutes Minimum	80 Minutes Minimum

¹Validated with KimGuard® KC600. (Validated to the following standards: ANSI/AAMI ST8:2008, ANSI/AAMI/ISO 17665-1:2006 and ANSI/AAMI ST79:2010.)

The reusable instruments should be visually inspected for wear and tear, routinely and prior to each surgical use. Visual inspection should be performed to look for imperfections or deterioration (such as corrosion, pitting, significant discoloration, cracking, breakage, disassembly, etc.). For reusable instruments that show evidence of such deterioration, the user should first contact Kinamed (or, for users located outside the USA, the authorized Kinamed distributor) to determine if repair is possible. If Kinamed determines that repair is possible, the user should thoroughly clean the deteriorated reusable instrument before returning for repair. If Kinamed determines that repair is not possible, then the deteriorated reusable instrument may be discarded using the hospital's normal disposal procedure for similar types of surgical instruments. Note that it is normal and expected for anodized instrument surfaces to show slight discoloration after repeated cleaning and sterilization cycles. Such discoloration is not evidence of deterioration.

CLEANING AND MAINTENANCE

All implants and instruments intended for end-user sterilization must be free of packaging material and biocleaners prior to sterilization. Cleaning, maintenance and mechanical inspection must be performed by hospital personnel trained in the general procedures of contaminant removal and use.

Pre-process Handling: Remove gross soiling by submerging the instrument into cold tap water immediately after use. Instruments should be cleaned within 30 minutes of use to minimize the potential of staining, damage, and drying.

(1) Manual Soak: Completely submerge instruments in neutral pH endozime detergent (i.e. Ruhof Endozime® AW Plus[®], dilute 1/2 ounce per one (1) gallon (4mL/liter) of warm (20-25°C) tap water) for 5 minutes. Use a soft bristled, nylon brush to gently scrub the device until all visible soil has been removed. Particular attention should be given to hard to clean areas. Visually inspect the instruments for cleanliness and repeat the manual soak step (1) until no visual contamination is present. For reusable instruments that are not visually clean after multiple attempts to repeat the manual cleaning cycle, the user may discard the instrument using the hospital's normal disposal procedure for similar types of surgical instruments. Note that it is normal and expected for anodized instrument surfaces to show slight discoloration after repeated cleaning and sterilization cycles. Such discoloration is not evidence of visual contamination.

Manual Rinse: Remove instruments from the enzymatic solution and rinse thoroughly under running Deionized (DI) water for a minimum of 1 minute. Thoroughly and aggressively brush and flush through cannulated areas using a DI water jet with the exit end submerged for a minimum of 1 minute. Flush instrument for a minimum of 2 minutes with DI water again to clear any debris.

(2) Combination Manual/Automated Soak and Rinse: Same as manual soak and rinse in section (1) above. **Washer/Dryer:** Place instruments in a suitable washer basket and load in an automatic washer/dryer. Cycle should be set for a Non Cautious wash cycle according to the table below.

Phase of Automatic Cycle:	Minimum Water Quality / Temperature	Minimum Time	Detergent / Dilution
Pre-Wash	Tap Water / 31°C	3 minutes	Not Applicable
Wash	Tap Water / 55°C	11 minutes	Ruhof Endozime® AW Plus [®] , pour directly into detergent reservoir at full strength or set the dispenser to 1/2 ounce per one (1) gallon (4mL/liter) of warm (20-25°C) tap water
Rinse	Tap Water / 34°C	1 minute	Not Applicable
Final Rinse -Thermal Disinfectant Rinse	DI Water / 94 °C	2 minutes	Not Applicable
Dry	Not Applicable / 121°C	15 minutes	Not Applicable
Total Cycle Time	Not Applicable	45 minutes	Not Applicable

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