KINAMED INCORPORATED	
Kinamed, Inc. 820 Flynn Road Camryn, CA 93012 USA 1-800-827-5775 www.kinamed.com	CE 0086 (FAX) 1-805-384-2792 ©2008 Kinamed™, Inc.
European Authorized Representative: MediMark® Europe Sarl. 11, rue Emile Zola B.P. 2332 F-38033 Grenoble Cedex 2, France Tel +33(0)476864322 Fax +33(0)476171982	B00109F
EN English SuperCable™ Iso-Elastic™ Cerglage System Instruction for Use	

INFORMATION CONCERNING USE OF THE SUPERCABLE™ ISO-ELASTIC™ CERCLAGE SYSTEM Protected by U.S. Patent Nos. 6,589,246 & 7,207,090. Other U.S. and International Patents Pending.

GENERAL REMARKS: This brochure is a general operative guide for implantation of the Kinamed SuperCable™ Iso-Elastic™ Cerglage polymer cable system. The Kinamed SuperCable™ Iso-Elastic™ Cerglage System consists of a braided cerclage cable and attached metal clasp. The cable is flexible and possesses high fatigue and tensile strength. The cable is made from biocompatible materials, consisting of UHMWPE strands braided over a nylon core. The clasp components are made from either stainless steel, titanium, or cobalt-chromium alloy (each package label states the material used for the enclosed clasp). The general principles of patient selection and sound surgical judgment apply to the cerclage procedure.

PACKAGING and LABELING: See Product Label for information regarding the specific product referenced in this document. The implant should be accepted only if received by the hospital or surgeon with the factory packaging and labeling intact.

CAUTION

This document sets forth recommended procedures for using Kinamed devices and instruments. It offers guidance, but, as with any technical guide, each surgeon must consider the particular needs of each patient and make appropriate adjustments as necessary.

CAUTION

Federal Law restricts this device to sale by or on the order of a licensed physician.

Ⓢ **Devices labeled for single use are intended to be used once only, for a single patient, because they may not perform as intended if they are reused. Reuse may lead to failure of the device to perform as intended.**

INDICATIONS: (1) Repair of long bone fractures due to trauma or reconstruction; (2) Reattachment of the greater trochanter of the hip with arthroplasty, surface replacement arthroplasty, or other procedures involving trochanteric osteotomy; (3) Stenotomy Closure; (4) Sublaminar and intrafascial wiring of the spinal column.

CONTRAINDICATIONS: (1) Active or suspected infection, either systemic or localized, in or around the implant site; (2) Patient conditions, mental or neurological, that would tend to impact the patient's ability to follow physician's instructions during the post-operative healing phase; (3) Skeletally immature patients; (4) Compromised vascularity that would inhibit adequate blood supply to the fracture or operative site; (5) Insufficient quality or quantity of bone, which would inhibit rigid device fixation; (6) Any disease affecting the support and fixation of the prosthesis; (7) Distant foci of infection, such as genitourinary, pulmonary, skin or other sites which may spread to the implant site. The foci of infection should be treated prior to, during, and after implantation.

STERILITY: Cable and clasp are supplied sterile. The package should be examined prior to use for possible breaks in the sterile barrier. Cable contains polyethylene and nylon polymers.

Do not autoclave or re-sterilize cable implants.

The SuperCable™ instruments are provided non-sterile and must be sterilized prior to surgical use per the following validated procedure:

Prevacuum/Pulsating:
 Wrapped, 3 – 3.5 minutes, 134°C – 137°C
 Sterilization times given above represent exposure time only and not total cycle time. (Validated to the following standards: FDA's 21CFR58, BS EN 554:1994, 556-1:2001 and ANSI/AAMI S779:2006)

PRECAUTIONS: The patient must be advised of the short and long term limitations of the procedure and the need to protect the implant from full weight bearing until adequate fixation and healing have occurred. Expected activity and inspection affecting the reattachment and/or fixation have been implicated in failure of this procedure by loosening, fracture, and/or wear of the implants. Loosening of the components can result in increased production of wear particles, as well as damage to the bone, making successful revision surgery more difficult. The patient should be cautioned to limit activities, protect the bone from unreasonable stresses, and to follow the instructions of the physician with respect to follow-up care and treatment. The patient should be warned of the surgical risks, and made aware of possible adverse effects. The patient should be warned that the device cannot and will not replicate the flexibility, strength, reliability or durability of normal healthy bone, that the implant can break or become damaged, particularly as a result of trauma, and that the device has a finite expected service life and may need to be replaced in the future. Transient bacteremia can occur in daily life. Dental manipulation, endoscopic examination, and other minor surgical procedures have been associated with transient bacteremia. To prevent infection at the implant site, it may be advisable to use antibiotic prophylaxis before and after such procedures. An implant should never be re-used. Even though an implant may appear undamaged, previous handling and in-service stresses may have created imperfections that would reduce the service life of the implant. Over-tensioning can create the potential for the cable to cut into the bone. Under-tensioning can cause inadequate fixation and increased fatigue stress. Improperly installed cable clasps can cause cable loosening. Excessively long cable ends caused by not cutting the cable flush against the clip may result in bursts. Proper handling of any implant is mandatory. Prior to surgical use, a visual inspection of each implant for possible imperfections should be performed. Damaged or altered implants may produce undesirable stresses and cause defects, which could lead to failure.

CAUTION

Cable and locking clasp may be re-tensioned and re-locked as necessary. If a locked cable must be removed from or relocated on the bone, it can be cut and discarded, and a new cable applied in the desired location. Alternatively, cable tension can be relieved by unseating the locking wedge with dedicated instrumentation. Intentional unseating of the wedge should be performed a maximum of once per cable.

WARNINGS: Improper placement, positioning and fixation of these devices may result in unusual stress conditions reducing the service life of the implant. Use care in utilizing the Cable Passer and other instrumentation to avoid damaging neurovascular structures. Exercise caution when positioning and tightening the Cerclage Cables to avoid impinging or tethering neurovascular structures and causing attendant significant sequelae. The surgeon is to be thoroughly familiar with the surgical procedure, instruments and implant characteristics prior to performing surgery. Periodic, long-term follow-up is recommended to monitor the position and state of the implants, as well as the condition of the adjoining bone. Patients should be questioned regarding sensitivity to metals. If there is a question pertaining to the patient's tolerance for titanium, stainless steel, or cobalt chrome, appropriate testing should be performed. The following conditions, singularly or concurrently, tend to impose severe loading on the affected bone, thereby placing the patient at high risk for failure: obesity, heavy labor, active sports participation, high general activity level, likelihood of falls, alcohol or drug addiction, and other disabilities. Some of the alloys used to produce orthopaedic implants contain metallic elements that may be carcinogenic in tissue cultures or intact organisms under unique circumstances. Questions have been raised in the scientific literature as to whether or not these alloys themselves may be carcinogenic in implant recipients. Studies conducted to evaluate this issue have not identified conclusive evidence of this phenomenon, in spite of the millions of implants in use. Polyethylene wear has been reported following implantation of many devices. Higher rates of wear may be initiated by metallic particles or other debris which can cause abrasion. Higher rates of wear likely will shorten the useful life of the device, and lead to an earlier than desired revision to replace the worn components. Cerclage cable breakage is reported in the literature. Possible causes of breakage include, but are not limited to, over tensioning the cable, fretting, and delayed or non-union of the fracture, osteotomy, or fusion site. In many instances, adverse effects may be clinically related rather than implant related. The surgical and post-operative management of the patient must be carried out with due consideration for all existing conditions of increased risk. Mental attitudes or personality disorders resulting in patient's failure to adhere to his/her physician's orders might delay post-operative recovery and aggravate adverse effects.

POST-OPERATIVE PROCEDURES: Late postoperative complications can include: (1) Early or late loosening, wear, and change in position of one or more components; (2) Trochanteric avulsion as a result of excess muscular tension, early weight bearing, or inadvertent intra-operative weakening; (3) Trochanteric nonunion due to inadequate reattachment and/or early weight bearing; (4) Progressive bone resorption and osteolysis. The patient should be cautioned to govern his/her activity level. Postoperative therapy should be structured to prevent excessive loading. The patient should be released from the hospital with complete instructions and warnings (preferably written) regarding further exercises and therapies. The patient should be encouraged to report any unusual changes to his/her physician.

RECOMMENDED CLEANING OF INSTRUMENTS: (Validated to NHS Health Technical Memoranda HTM2030 – Washers and HTM2010 – Sterilizers)
(1) Manual Soak: Completely submerge instruments in neutral pH endozime detergent for 5 minutes. Use a soft bristled, nylon brush to gently scrub the device until all visible soil has been removed. Particular attention should be given to hard to clean areas.
Manual Rinse: Remove instruments from the enzymatic solution and rinse thoroughly under running tap water. Thoroughly and aggressively brush and flush through canulated areas using a water jet with the exit end submerged. On the Tensioning Instrument turn handle counter-clockwise to fully open the cable gripping jaws. Flush instrument, including jaws, to clear any debris. Next, turn knob of tensioning instrument to fully expose remainder of threads on the lead screw. Flush instrument again to clear any debris.

(2) Combination Manual/Automated Soak and Rinse: Same as manual soak and rinse in section (1) above. **Washer/Drier:** Place instruments in a suitable washer basket and load in an automatic washer/drier. Cycle should be set for a Non-Cautish wash cycle for a duration of 70 minutes using a neutral pH endozime detergent. The endozime detergent should be used at a specified concentration in a 14-minute cleaning cycle.
MAINTENANCE OF INSTRUMENTS: Prior to autoclave sterilization, apply a surgical grade lubricant to the lubricant and the mechanical surfaces of the tensioning instrument. Be sure that the lubricant fully penetrates the mechanism. Wipe any excess lubricant that may have been deposited on the back of the instrument body. Extreme care should be

exercised lubricating this thread to prevent any excess lubricant from depositing in the cable holding grooves.

CAUTION

Use extreme care to prevent excess lubricant from reaching the cable holding grooves. Lubricant in cable holding grooves may cause cable slippage during tensioning.

CARE and HANDLING: Use extreme care in handling and storage of implant components. Cable and clasp must be handled with care. Twisting, kinking, cutting, notching, or scratching the braided surface may reduce the strength, fatigue resistance and/or wear characteristics of the implant system. These, in turn may induce internal stresses that are not obvious to the eye and may lead to fracture of the component. Implants and instruments should be protected during storage from corrosive environments, such as salt air, etc. Only instruments designed for use with this system should be used to assure correct implantation. Review of these handling instructions is important. Damaged instruments may lead to improper cable tension or implant position and result in implant failure. Thorough familiarity with the surgical technique is essential to ascertain their proper working condition.

ADVERSE EFFECTS: With all implants, asymptomatic localized progressive bone resorption (osteolysis) may occur around or remote from the prosthetic components as a consequence of foreign body reactions to particulate matter. Particulate matter is generated by the interaction between implant components, as well as between the component and bone, primarily through wear mechanisms of adhesion, abrasion, and fatigue. Osteolysis can lead to future complications necessitating the removal and replacement of prosthetic components. Although rare, metal sensitivity reactions in patients following device implantation have been reported. Implantation of foreign material in tissues can result in cellular reactions involving lymphocytes, macrophages and fibroblasts. Early, or late loosening, cracking, fracture or deformation of one or more components. This is often attributable to factors listed under "WARNINGS." Loosening may also occur due to defective fixation or improper positioning. Early or late infection, which may require removal of the implant.

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