



NeuroPro[®]

Chiari Plate

Pre-Shaped Mesh Design

Decompression Protection




Effective and Affordable

The Chiari Plate (part no. 14-1P3-0700, illustrated below) is used to cover the decompressed posterior fossa in the treatment of Chiari malformation. It provides an efficient means to cover this area with a titanium mesh to maintain the decompression and to prevent adhesions of muscle to dura which may cause postoperative patient headaches. This Chiari plate is pre-shaped to cover the typical decompressed area.



Implant

Catalog #	DESCRIPTION
14-1P3-0700	Chiari Plate



KINAMED
INCORPORATED

For more information:

Phone (805) 384-2748
 Toll-Free (800) 827-5775
 Fax (805) 384-2792
 Website www.kinamed.com

820 Flynn Road, Camarillo, CA 93012-8701



ISO 13485
FM 75124

