

# VacuJet®



Instantly clear clogged suction

A better way to clean fluids and debris

Does this look familiar?




A rapid blast of CO<sub>2</sub> from the VacuJet quickly and efficiently clears a clogged suction tip



Self-Cleaning Suction

Save valuable O.R. Time<sup>1</sup>

Eliminate mess and cost



"For years, I have wanted a simple elegant suction tool that eliminated the fiddle of chronic plugging. Kinamed has created the perfect solution with a sucker on the pistol grip CarboJet handpiece."

**James B. Stiehl, M.D.**  
Centralia, Illinois, USA



**KINAMED**  
INCORPORATED



*Expect Innovation.*

# VacuJet® Self-Cleaning Suction Enhancement

- Traditional suction wands are easily clogged
- Cleaning/Unclogging wastes valuable O.R. time and increases procedure cost
- Clogged suction is a distraction for surgeon & staff
- The VacuJet is a nozzle “add-on” for the CarboJet® system that includes a suction tube with a filter screen at its distal end
- A rapid blast of CO<sub>2</sub> quickly and efficiently clears the blockage allowing the surgeon to get back to work with minimal disruption

**In a prospective randomized study, VacuJet was shown to save significant O.R. time versus two different traditional suction devices<sup>1</sup>**

## VacuJet® Items

<b>Catalog</b>	<b>Description</b>	
25-100-2001	CarboJet CO <sub>2</sub> Tube Set, Sterile Packed (each)	
25-200-0200	CarboJet Handpiece	
25-500-0242	VacuJet Suction Tube	
25-500-0244	VacuJet CO <sub>2</sub> Nozzle	
25-500-0246	VacuJet Suction Tube, Long	
25-500-0248	VacuJet CO <sub>2</sub> Nozzle, Long	

1. Stiehl (2014) Mechanical Performance of a Self-Unplugging Surgical Suction Instrument: A Randomized Controlled Trial. *Reconstructive Review* 4(1):18-22.



**K I N A M E D**  
INCORPORATED

*Expect Innovation.*

**For more information:**

**Phone** (805) 384-2748  
**Toll-Free** (800) 827-5775  
**Fax** (805) 384-2792  
**Website** [www.kinamed.com](http://www.kinamed.com)

**820 Flynn Road, Camarillo, CA 93012-8701**



ISO 13485  
FM 75124



To view a video demonstration, visit us on the Web at: [www.kinamed.com](http://www.kinamed.com)